
Bedienungsanleitung · Operating instructions · Istruzioni d'uso
Mode d'emploi · Gebruiksaanwijzing · Manual de instrucciones

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Digital Incubator for Reptile Eggs

Thank you for purchasing the Easy Breeder. For ideal hatching rates and healthy young animals, please read these operating instructions before initial operation.

Equipment: The Easy Breeder is supplied with a two-part housing of polystyrene. This is only designed for safe transport of the eggs and is not needed for normal operation. The incubator housing itself consists of a transparent bottom section of plastic, allowing visual inspection of the eggs in the interior from all sides. The digital controls, the display indicating the current temperature, relative humidity and a day counter showing how many days your eggs have been in the Easy Breeder are integrated in the cover. The modern combined heating-ventilation system is a technical highlight. This allows incubating the eggs with temperatures hardly fluctuating and ensures a first class hatching rate! The sensors for relative humidity and temperature are safely accommodated under a grid in the cover. Water channels are provided in the floor for regulating the humidity during the breeding and hatching phases. The second grid included in the delivery is used for incubation on two levels, so that twice the number of eggs can be incubated.

Table of contents:

1. Starting up:

2. Setting the temperature

3. Alarm setting for temperature
4. Alarm setting for relative humidity
5. Defining a restricted temperature range
6. Resetting the day counter / Reset to factory setting
7. Useful hints/Incubation table
8. Cleaning and care of the Easy Breeder
9. Problems and solutions
10. Safety instructions and technical data
11. Declaration of conformity
12. Warranty

1. Starting up:



- 1.1 First connect the power cable to the Easy Breeder and then plug the connector into the mains socket.
- 1.2 Now press the  **ON / OFF - key** to switch on the Easy Breeder
- 1.3 An alarm signal sounds. **You can clear alarms at any time by pressing the   arrow keys.**
- 1.4 Open the Easy Breeder fill the 7 water channels using the included bottle. These are located underneath the deposition grid. Each channel has a capacity for 1 full bottle (approx. 200 ml) of water.
- 1.5 Set the temperature to the desired value. Different temperatures are required for different reptiles (see incubation table).

2. Setting the temperature:



- 2.1 Press the  **menu - key** control panel once.
- 2.2 Press the  **up arrow -** or  **down arrow - key** to set the desired temperature.
- 2.3 Press the  **menu - key** again to complete this process.
The temperature indicator in the control panel always shows the current temperature. The  **status light** will show when the Easy Breeder is activated. The setting range is 20.0 – 39.5° C. Adjustable in 0.1° C increments.

3. Alarm setting for temperature: The factory setting for the temperature alarm is programmed at 1° C above or below the defined desired temperature. When this value is exceeded an alarm signal sounded. If you wish to change this setting, proceed as follows:

3.1 Setting the low temperature alarm (AL):

- 3.1.1 Hold the  **menu - key** for approx. 3 seconds.
- 3.1.2 Press the  **up arrow -** or  **down arrow - key** until the display reads (AL).
- 3.1.3 Press the  **menu - key** again to confirm.
- 3.1.4 Use the   **arrow keys** to set the lower alarm limit. The setting range is 0.0 – 15.0° C. Can be set in 0.1° C steps.
- 3.1.5 Press the  **menu - key** again to complete the process.

3.2 Setting the high temperature alarm (AH):



- 3.2.1 Hold the  **menu - key** for approx. 3 seconds.
- 3.2.2 Press the  **up arrow -** or  **down arrow - key** until the display shows (AH).
- 3.2.3 Press the  **menu - key** again to confirm.
- 3.2.4 Use the   **arrow keys** to set the lower alarm limit. The setting range is 0.0 – 15.0° C. Can be set in 0.1° C steps.
- 3.2.5 Press the  **menu - key** again to complete the process.

4. Alarm setting for relative humidity (AS): The factory setting for the  relative humidity alarm is set to 45 %. When the humidity falls below this value, an alarm signal sounded. If you wish to change this setting, proceed as follows:

- 4.1 Hold the  **menu - key** for approx. 3 seconds.
- 4.2 Press the  **up arrow -** or  **down arrow - key** until the display shows (AS).
- 4.3 Press the  **menu - key** again to confirm.
- 4.4 Use the   **arrow keys** to set the desired alarm limit. The setting range is 1 – 99 %  relative humidity. Can be set in steps of 1 %.
- 4.5 Press the  **menu - key** again to complete the process.

5. Defining a restricted temperature range: This function prevents that the temperature accidentally exceeds or drops below the defined desired temperature. The factory setting has no defined limit. If you wish to change this setting, proceed as follows:

5.1 Setting the upper temperature limit (HS):

- 5.1.1 Hold the  **menu - key** for approx. 3 seconds.
- 5.1.2 Press the  **up arrow -** or  **down arrow - key** until the display shows (HS).
- 5.1.3 Press the  **menu - key** again to confirm.
- 5.1.4 Use the   **arrow keys** to set the desired upper limit. The setting range is 20 - 50° C. Can be set in 0.1° C steps.
- 5.1.5 Press the  **menu - key** again to complete the process.

5.2 Setting the lower temperature limit (LS):

- 5.2.1 Hold the  **menu - key** for approx. 3 seconds.

5.2.2 Press the  **up arrow** - or  **down arrow** - key until the display shows (LS).

5.2.3 Press the  **menu** - key again to confirm.

5.2.4 Use the   **arrow keys** to set the desired lower temperature limit. The setting range is 20 – 50° C. Can be set in 0.1° C steps.

5.2.5 Press the  **menu** - key again to complete the process.

 **For example:** If you have defined an upper temperature limit of 28.2° C and a lower limit of 27.4° C, you will only be able to set the temperature – as described in item 2 – within this range.

6. Resetting the day counter:

6.1 The Easy Breeder has an integrated day counter (called  **“Day”** in the display). This means you can always see how long the eggs have been in the incubator. To reset the day counter, press and hold the  **Reset** - key a few seconds. This will cause the counter to jump back to **“0”**.

7. Useful hints / Correct handling of the eggs / incubation table

7.1 Useful hints:

- Operate the Easy Breeder in a room with a constant temperature, because external temperature fluctuations can affect the incubator. It should be set up away from heat sources and direct sunlight. Make sure it rests on a horizontal base and is free from vibrations or shock and is free-standing without obstructions to the air flow.
- It may take several hours until the temperature stabilises at the desired level, depending on room and set temperature, sensor position and the number of bowls with incubation material. We recommend you operate the incubator for one day without hatching eggs so that the incubation material can fully absorb the heat.
- After this day of operation for the bowls with incubation material, add the eggs. Opening the Easy Breeder and the number as well as the intrinsic temperature of the hatching eggs can cause a temperature drop, a gradual increase and temperature fluctuations. Please do not change the programming during this period!
- During the incubation period, minor temperature and humidity fluctuations may also occur. This is caused by the room temperature, the number of incubation bowls and eggs, the incubation material and the relative humidity in the room. This has no harmful effect on the incubation process. Constant temperatures do not occur in natural environments either.
- In principle, the Easy Breeder should only be opened when needed (to inspect the eggs or replenish liquid) to ensure a constant climate.

 **Expert hint:** Even when incubating only one clutch of eggs the incubator should be equipped with cricket boxes or other containers and filled with moist substrate such as vermiculite , soil or similar material. This significantly increases the heat capacity of the entire system keeps temperature fluctuations low. The incubator will also maintain its temperature much longer in the event of a power failure.

7.2 Correct handling of the eggs:

- Incubating the eggs in the Easy Breeder has the purpose of achieving as high a hatching rate as possible under controlled and optimum conditions and protecting the eggs against external influences.
- Only collect the eggs once the female has completely finished depositing the eggs. This will ensure that the female is not subjected to stress, which may cause a fatal retention of eggs.
- After egg deposition, remove the eggs from the terrarium. It is essential that the eggs are placed into the incubation substrate (e.g.: Vermiculit item no.: 36320 and 36325) in the same position as before their removal to prevent the poles of the eggs from shifting (reptile eggs do not have any retaining chords („chalazas“) like birds have). This could cause the foetus to die.
- If you wish to mark the eggs, wait until they have fully hardened to avoid damaging the egg membrane. All eggs should always be handled with care.
- The moisture of the incubation substrate plays a key role in the development of the reptile eggs. Here the rule applies that soft-shelled eggs have a higher moisture requirement than hard-shelled ones. Too little substrate moisture would extract moisture from the eggs and this loss of moisture would be fatal for the developing animal.
- We can recommend vermiculite (e.g.: Vermiculit item no.: 36320 and 36325) as a good incubation substrate. With this substrate it is possible to regulate the substrate moisture by means of adding water and the grain size. fine grained vermiculite (e.g.: Vermiculit 0 – 4 mm item no.: 36320) is recommended for small and soft-shelled eggs. Coarser vermiculite (e.g.: Vermiculit 3 – 6 mm item no.: 36325) is recommended for large and hard-shelled eggs.
- Regularly check the moisture of the incubation substrate. When required, dose the water directly into the incubation substrate and not above the eggs.
- The time of hatching can differ greatly, even within the same clutch of eggs. Primarily this is caused by differences in temperature and humidity.

7.3 Incubation table: The table below contains the guidance values for a number of reptile species based on experience. These do not guarantee successful breeding!

		Incubation temperature in °C	Duration of incubation:
Common agama	Agama Agama	30	approx. 60 – 70 days
Green ameiva	Ameiva Ameiva	30	approx. 130 days
Red-throated anole	Anolis carolinensis	27 – 30	approx. 40 – 60 days
Knight anole	Anolis equestris	27 – 30	approx. 70 – 100 days
Bahaman anole	Anolis sagrei	27 – 30	approx. 40 – 60 days
Plumed basilisk	Basiliscus plumifrons	30	approx. 55 – 65 days
Striped basilisk	Basiliscus vittatus	30	approx. 50 – 70 days
Fischer's chameleon	Kinyongia fisheri	increasing from 17, 22 and 25	approx. 300 days
Killimanjaro two-horned chameleon	Kinyongia tavetanum	day 22 – 23, night 18	approx. 140 days

Veiled chameleon	Chamaeleo c. calypttratus	28	approx. 180 days
Helmeted iguana	Corytophanes cristatus	28	approx. 100 days
Desert collared lizard	Crotaphytus bicinctores	28 – 32	approx. 45 – 60 days
Eastern collared lizard	Crotaphytus collaris	28 – 32	approx. 45 – 60 days
Black iguana	Ctenosaura similis	28 – 29	approx. 80 – 100 days
Pegu forest gecko	Cyrtodactylus peguensis w	25 – 28	approx. 90 – 95 days
African egg-eating snake	Dasyplectis scabra	30	approx. 60 – 70 days
Desert iguana	Dipsosaurus dorsalis	28 – 33	approx. 45 – 80 days
Russian ratsnake	Elaphe schrenckii	28	approx. 40 – 60 days
Leopard gecko	Eublepharis macularius	28	approx. 45 – 55 days
Panther chameleon	Furcifer pardalis	26 – 28	150 – 300 days
Takay	Gekko gekko	30	approx. 100 days
Lined gecko	Gekko vittatus	28 – 30	approx. 90 – 160 days
Chinese leopard gecko	Goniurosaurus luii	26 – 28	approx. 70 – 80 days
African fat-tailed gecko	Hemitheconyx caudicinctus	28	approx. 70 – 90 days
Green iguana	Iguana iguana	29	approx. 70 – 90 days
Californian kingsnake	Lampropeltis g. californica	28	approx. 60 – 70 days
Pueblan milk snake	Lampropeltis t. campbelli	28	approx. 55 – 65 days
Northern curly-tailed lizard	Leiocephalus carinatus	28	approx. 60 – 70 days
Hispaniolan masked curly-tailed lizard	Leiocephalus personatus	28	approx. 60 – 70 days
Reeve's butterfly lizard	Leiolepis reevesii	28	approx. 60 – 70 days
White-lipped python	Leiopython albertisii	29 – 31	approx. 55 – 70 days
Carpet python	Morelia spilota variegata	29	approx. 70 – 85 days
Green tree python	Morelia viridis	29 – 30	approx. 45 – 55 days
Rough green snake	Ophedrys aestivus	27 – 28	approx. 50 – 60 days
Great plains grass snake	Pantherophis emoryi	27	approx. 60 days
Corn snake	Pantherophis guttatus	27	approx. 60 days
Desert horned lizard	Phrynosoma mclartyrhinos	27 – 28	approx. 45 – 55 days
Chinese water dragon	Physignathus cocincinus	28	approx. 70 – 90 days
Rankin's dragon	Pogona henrylawsoni	30	approx. 60 – 70 days
Central bearded dragon	Pogona vitticeps	28,5 – 30	approx. 60 – 70 days
Blood python	Python curtus	29 – 31	approx. 70 days
Burmese python	Python molurus bivittatus	29 – 31	approx. 60 – 90 days
Royal python	Python regius	29 – 31	approx. 70 days
Crested gecko	Rhacodactylus ciliatus	36 – 28	approx. 60 – 120 days
Chuckwalla	Sauromalus obesus	31	approx. 70 days
Rosebelly lizard	Sceloporus variabilis	28	approx. 50 – 60 days
Six-striped long-tailed lizard	Takydromus sexlineatus	24 – 25	approx. 30 – 45 days
Viper gecko	Teratolepis fasciata	26	approx. 60 – 70 days
Dwarf wonder gecko	Teratoscincus microlepis	28	approx. 85 – 150 days
Red eyed crocodile skink	Tribolonotus gracilis	24	approx. 60 – 70 days
Four horned chameleon	Triceros q. quadricornis	20 – 25	approx. 125 – 160 days
Argentine black and white tegu	Tupinambis merianae	30	approx. 70 – 80 days
Mop head iguana	Uranoscodon superciliosus	26	approx. 90 – 110 days
Ocellated uromastyx	Uromastyx ocellata	32	approx. 80 – 90 days
Ornate uromastyx	Uromastyx ornata	32	approx. 80 – 90 days
Side-blotched lizard	Uta stansburiana	28	approx. 40 – 45 days
Savannah monitor lizard	Varanus exanthematicus	27 – 31	approx. 115 – 180 days
Shield tailed agama	Xenagama taylora	29	approx. 50 – 55 days

For more detailed information we recommend the book: "Die Inkubation von Reptilieneiern" (The Incubation of Reptile Eggs by Dr. Gunther Köhler published by Herpeton Publishing.

8. Cleaning and care of the Easy Breeder: The Easy Breeder should always be thoroughly cleaned using a damp cloth and and disinfected with Rextix Desinfect (item no.: 38002) after each incubation.

9. Problems and solutions: When you encounter problems, check you settings on the Easy Breeder and carry out a reset if necessary (disconnect the mains plug and plug back in while pressing and holding the **Reset** button). Should this not lead to success, please speak to your specialist dealer.

10. Safety instructions and technical data

10.1 Safety instructions:

- The product must only be operated as supplied
- The device must not be operated if there are any defects on feed lines, heating element or display.
- Electronic components must not be modified in any way, and the leads must not be shortened.
- Disconnect the mains plug prior to opening, maintenance or cleaning.
- The heating element must be cleaned of dust and dirt particles before initial use.

- The heating element must only be touched after it has cooled down to avoid burns.
- The top section with the electrical installations must only be cleaned with a damp cloth.
- Do not immerse in water or hold into a jet of water.

10.2 Technical data:

Temperature display range:	-9.9° C – 99.9° C
Temperature control range:	20.0° C – 39.5° C
Temperature control accuracy:	0,1° C
Display accuracy:	0,1° C
Programming and measurement range for alarm:	0 – 15° C
Heating element output:	80 W
Operating voltage:	240 V ~

11. Declaration of conformity: The manufacturer certifies that the devices conform with EMC Directive 2014/30/EU, LVD Directive 2014/35/EU & ROHS Directive 2011/65/EU.

12. Warranty: You can find the warranty card with the warranty conditions in the Service section of our website: www.dohse-terrarium.com

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Dohse Aquaristik GmbH & Co. KG